

# Instrument disposal for the medi-matic® 115.7 series examination and treatment chair for gynaecology



**New**

## **Instrument disposal Model No. 101.204**

Electromotive, of sheet steel, powder coated in the colour of the chair, fitted on the column, connection cable with mains plug, connection: 90 – 240 V, 50/60 Hz, power consumption max. 20 W

Illustration:  
medi-matic® 115.755 with instrument disposal 101.204, foot supports 101.442, leg plate 101.445 and head cushion 101.372  
Powder coated finish grey aluminium RAL 9007, upholstery gold yellow 86, plastic parts silver grey

# Detailed functionality



Electromotive opening and closing by pressing the front panel for instrument disposal, fully extendable on telescopic rails, frame with removable 12 litre stainless steel basin (for 5 litres of liquid), plastic screen insert to collect about 15 specula. Additional stainless steel storage area. Maximum load 8 kg. Can be retrofitted (series 115.7) on request.

## Note:

Power supply via separate mains connection (cables of 2 m + 5 m).

## Getting to the point – Your key benefits:

- Instrument disposal with direct access
- Easy to use
- Space-saving solution
- No additional disposal module required
- Matching in colour and design

## Can also be used in conjunction with:

Colposcope support 101.430/101.435 **or**  
Schmitz Set with supporting arm 101.650  
**or** shelf with supporting arm 101.651

## Cannot be used in conjunction with:

Footrest 101.475  
Mounting step 101.452  
Stainless steel basin 101.201/101.202



Opening and closing by pressing



Electromotive opening and closing



Removable stainless steel basin with screen insert