

varimed® Height adjustable Examination and Massage couches




NEW




- Large and ergonomic height adjustment range
- No displacement of table-top during height adjustment
- Moulded base provides generous legspace



varimed Examination and Massage couch €



211.410



211.420 with 200.025



1



2



3

Description

Examination and Massage couch € 211.410

hydraulic height adjustment 500 to 920 mm by means of foot lever. Servo-assisted head section adjustment from -20° to +40°. 2-section upholstery 70mm thick, leatherette cover acc. to Schmitz colour chart. Powder coated chassis in white aluminium RAL 9006, grey aluminium RAL 9007 or grey white RAL 9002. Silver-grey plastic chassis cover. With 4 levelling feet, not electrically conductive. Safe working load: 225 kg

Examination and Massage couch € 211.420

as 211.410, but electromotive height adjustment from 500 to 920 mm by means of foot control unit, with safety ON/OFF switch to turn control element on or off, 2.40 m connecting cable with mains plug. Upon request, hand control instead of foot control unit.

Accessories:

1

Mobility 200.025

4 closed double castors Ø 100 mm, 2 of which el. conductive, central lock per axis, 1 castor with directional lock (Height adjustment range remains the same)

2

Side guards, pair 200.040

length: 990 mm, folding down to below top of mattress, powder coated white aluminium RAL 9006

3

Holder for foot control 200.035

for Examination and Treatment couch 211.420, of steel, powder coated white aluminium RAL 9006, for mounting at the long side of the couch (retrofitable)

Paper roll holder 200.093

for paper roll up to 590 mm wide

Paper roll 200.103

500 mm wide

Edge protection 200.041

of plastic, grey, at all four upholstery corners. (w/o illustration)

Head section with nose outlet 200.050

Surcharge for coloured chassis 200.653

acc. to colour chart

*** 4 Colours for leatherette upholstery in "classic" quality without surcharge
13 Colours for leatherette upholstery in "comfort" quality (surcharge 200.852).
Please see pages 26 and 27, catalogue no. 70.**

**Standard powder-coating colours for frame: RAL 9002, RAL 9006 or RAL 9007,
Please indicate desired colour upon order placement.**

The company reserves the right to alter the design and measurements without notice.

Schmitz Catalogue No. 19, 11.2011

Schmitz u. Söhne GmbH & Co. KG, Zum Ostenfeld 29, 58739 Wickede (Ruhr), Germany
Telefon +49 (0) 2377 84 0, Telefax +49 (0) 2377 84 162, www.schmitz-soehne.com