

PRESS-RELEASE

More Convenience for the Gynaecological Workplace

Space-saving Instrument Container for the medi-matic®-Series

Wickede/Duesseldorf, November 2009. The new instrument container for the medi-matic®-series is a remarkable labor-saving device for gynaecologists. The container can be operated with great ease and convenience. Its space-saving integration and easy access make an additional disposal unit for the gynaecological workplace superfluous. Moreover, the container can be adapted in colour and design to the workplace and the gynaecological chair. By touching the faceplate the container opens and closes electromotively. It consists of a frame with a removable 12 litre SST container, a plastic filter insert for disposal of 15 specula and an additional storing tray also made of SST. The container can be fitted subsequently to all models of the medi-matic 115.7 series.

hall 13 stand 13 A 56