

# Transfer mattress for easy and comfortable transfer of patients for patient stretchers STS 282 and STL 285



80 mm thick, of resistant, hard-wearing PU material, antistatic, black, suitable for X-ray imaging, underside of high quality special sliding material, 1,950 x 600 mm, with 4 looped handles along each side. To be secured on both sides of the stretcher by Velcro straps.



**Transfer mattress, Model No. 206.389**  
for patient stretcher STS 282.101

**Transfer mattress, Model No. 207.389**  
for patient stretcher STL 285.080 and 285.100

# Getting to the point – Your key benefits



Secured easily and safely by Velcro straps on both sides



- High sliding quality and optimum comfort
- Less strenuous and easy to use
- Helps to save time and costs as it can be handled by just 1 person
- Ideal for transfer to X-ray tables
- Suitable for patients up to 225 kg
- No additional investment, as it can also be used as a standard mattress

## **Important:**

The transfer mattress may only be used on Schmitz patient stretchers 282.101, 285.080 and 285.100 manufactured after 02/2009!

The transfer mattress cannot be secured on previous models nor is it suitable for use on patient stretchers from other manufacturers!